

About InnAccel

About InnAccel – InnAccel is a Bangalore-based company that has created India's first MedTech innovation platform, anchored by state-of-the-art R&D facilities and a network of national and global partners. Our goal is to build a portfolio of medical products designed, engineered, and priced, for global emerging markets.

 +91 80 40923864

 info@innaccel.com

 www.noxeno.in

Made possible by the generous support of:



सत्यमेव जयते

DEPARTMENT OF BIOTECHNOLOGY
GOVERNMENT OF INDIA

Ministry of Science and Technology, Government of India

SCHOOL OF INTERNATIONAL
BIODESIGN



nōxeno

A Simple and Safe Nasal Foreign Body Extractor



InnAccel
innovation... accelerated

BENEFITS

Easy

- Simple, two-step assembly.
- Ergonomic, single-handed use for fast NFB removal.



Safe

- Atraumatic tip design with soft sheath.
- Integrated, tip-focused light for visualization.



Durable


- Modular and replaceable components.
- Robust component design and materials.

PRODUCT FEATURES



Handle

- 5mm LED light
- Light button
- Replaceable batteries (2x AAA)



Extractor Assembly

- Optimal tip size and angle
- All Stainless-steel
- Strong hinge for pull forces up to 20N
- Modular for easy reprocessing
- Autoclavable



Rubber Sheath

- Single-patient use disposable
- Sterile
- Soft for patient comfort